CellaLab
Laboratory Information
Management System
für Labordatenanalyse

CellaLab
Laboratory Information
Management System
Information

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Introduction

CellaLab was developed as a data management solution for laboratories to monitor samples and support product quality control. For more than 30 years, KELLER HCW GmbH has been providing control systems for the chemical industry and food technology sector. Based on this experience and expertise, KELLER has designed yet another enterprise software application, employing technologies which have proven highly robust and reliable and extending its wide range of process control solutions. CellaBatch features client/server architecture based on Microsoft .NET-Framework 3.5. The features

of the integrated SQL Server 2008 create a wide variety of capabilities, utilizing advanced bi-functionality as part of a comprehensive bi-platform.

CellaBatch can be used as a stand-alone solution or integrally with KELLER's CellaTrace process control system. When used together with CellaTrace, you will not need to define the products and the process cells in CellaLab, since this data will be communicated from a mutual master data base, thus preventing double entries or inconsistencies.

Stamm Master Stammdaten data


Master Data

All data relevant to hazard analysis and control point monitoring is entered into the master data and managed there. Information such as the definitions of article, client and analysis are part of this data. The real material analyses communicate with this master data. If a connection exists with CellaTrace, then the "article" will be imported from the CellaTrace data base.

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
The sample test results are stored in “Analysis definition” and available for further processing. You can add various details to the stored “Analysis definition” such as the defined measurand, the reference value (optimum) and, if required, a specific calculation formula. By factoring in the reference value, the analysis results are checked against the upper and lower threshold values of the applied profile. The terms used for the individual methods as well as the terms for the individual measurands are stored in charts which can be accessed and augmented by the user.

With CellaLab, the employed laboratory apparatus as well as the containers involved in sample testing can be defined. Several containers can be assigned to one apparatus. Each container, in turn, can be assigned to a specific analysis result. Therefore it is possible to manage the sample data for various samples which have been tested by one apparatus. Barcode scanner integration is an option.

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CellaLab enables the user to create sample types. Each sample type can have analysis definitions assigned to it. The user can create and custom define every type of sample.

CellaLab automatically generates the analyses to be performed when a new material sample of a specific sample type is created.

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
A standard article profile can be created in CellaLab's Profile Maintenance. In the respective testing scheme the appropriate definitions and their respective threshold values can be established. If desired, an optimum value can be entered. The permissible attribute value can also be entered for analyses denoted by a clear yes/no finding.

Depending on the article profile, any number of user-specific profiles can be created. The testing scheme can be copied from other product profiles and adapted or augmented as required. If the limit value – which has been defined for the test criteria of the superordinate standard product profile – is not defined within the client-specific profile, then CellaLab will check the limit value data against the superordinate product profile.

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
Data Collection


Before test data is col-
lected, a material data set
is created. This data set
contains information about
the article and particulars
on the material. The mate-
rial data set automatically
draws on the standard
product profile as a basis
for comparison. Any num-
ber of samples can be
assigned to a material data
set. The test data for each
of these analyses are
entered here.

Before creating a sample,
the material must be
selected to which the
sample will be assigned.
The user can either create
a new material data set or
use an existing data set.

Analysen Analyses of material samples

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For each material sample, any number of analyses can be stored. When a new sample has been created, and based on the selected type, those test results which apply to this sample type will be generated. The status of the individual analysis and any possible limit violations (based on the profile) will be shown.

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

All yielded analyses for a particular material are summarized in the evaluation. In order to generate this evaluation, first the specific material must be selected. The analyses for all material samples are all available to choose from. This selected sample data is used for the evaluation. If several test results are selected for a specific analysis definition, CellaLab will calculate an average value. In case a profile has been assigned to the material, the results of the evaluation are checked against the respective testing scheme and shown highlighted. The results can be displayed/printed in report format.

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Three companies – one strategy: Professionals in heavy clay works

KELLER HCW – for more than 100 years one of the worldwide leading machine and plant constructors.

Starting with a cutter, followed by dryers and kilns to packaging plants, the complete handling, automation and Measuring & Controlling.

The **Rieter** works – for more than 100 years specialized in machines and plants for preparation and shaping – from storage systems to feeders, crushers, mixers, batching equipment to extruders.

novoceric – one of the leading suppliers of grinding machines for calibrated bricks with long lasting experiences in the field of handling and handling plants.

As professionals in heavy clay works these three companies are your qualified partners – with innovative new plants as well as with the modernization and optimization of existing plants.

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KELLER HCW GmbH · Carl-Keller-Straße 2-10 · 49479 Ibbenbüren-Laggenbeck · Germany · Telefon: +49 54 51 85-0